

STAR Protocols

Building a “Hello World” for Self-driving Labs: The Closed-loop Spectroscopy Lab Light-mixing Demo (CLSLab:Light)

--Manuscript Draft--

| Manuscript Number: | STAR-PROTOCOLS-D-23-00017R3 |
|---|---|
| Full Title: | Building a “Hello World” for Self-driving Labs: The Closed-loop Spectroscopy Lab Light-mixing Demo (CLSLab:Light) |
| Article Type: | Protocol |
| Corresponding Author: | Sterling Baird The University of Utah Salt Lake City, UT UNITED STATES |
| Corresponding Author Secondary Information: | |
| Corresponding Author's Institution: | The University of Utah |
| First Author: | Sterling G. Baird |
| Order of Authors: | Sterling G. Baird Taylor D. Sparks |
| Abstract: | Learn how to build a Closed-loop Spectroscopy Lab: Light-mixing Demo (CLSLab:Light) to perform color-matching via RGB LEDs and a light sensor for under 100 USD and less than an hour of setup. Our tutorial covers ordering parts, verifying prerequisites, software setup, sensor mounting, testing, and an optimization algorithm comparison tutorial. We use secure IoT-style communication via MQTT, MicroPython firmware on a pre-soldered Pico W microcontroller, and the self-driving-lab-demo Python package. A video tutorial is available at https://youtu.be/D54yfxRSY6s . |
| Additional Information: | |
| Question | Response |
| Original Code<p>Does this manuscript report original code? | Yes |
| Reviewers must have anonymous access to these original code that is free-of-cost. Please provide code location and instructions for access here.<p>Please consult this Author's guide for more information: How standardized datasets and original code accompany Cell Press manuscripts from submission through publication" or email us at joule@cell.com. as follow-up to "Original Code<p>Does this manuscript report original code?" | https://github.com/sparks-baird/self-driving-lab-demo |
| Standardized datasets<p>A list of datatypes considered standardized under Cell Press policy is available here. Does this | No |

manuscript report new standardized
datasets?



Materials Science and Engineering

122 S. Central Campus Drive, Salt Lake City, Utah 84112 (801) 581-8632

December 20, 2022

Dear Dr. Leila Shokri:

We are very pleased to submit this invited protocol *Build instructions for Closed-loop Spectroscopy Lab: Light-mixing Demo*, to Star Protocols as a companion manuscript to our Matter perspective article: *What is a minimal working example for a self-driving laboratory?* This work presents the build instructions for a self-driving laboratory demonstration that costs under \$100 and requires minimal setup space and time.

The demo is a minimal working example (MWE) for a materials acceleration platform that meets the three requirements of a MWE—minimal, complete, and reproducible—and is suitable for low-cost prototyping of materials acceleration platform (MAP) concepts such as distributed autonomous networks and as an educational tool to increase hands-on exposure to setting up MAPs.

Additional software, documentation, and build instructions are being developed on GitHub and Hackaday. This manuscript represents a fixed and reproducible snapshot.

Sincerely,

A handwritten signature in black ink, appearing to read "Taylor Sparks".

Dr. Taylor Sparks
Associate Professor & Associate Chair
Materials Science and Engineering Department
University of Utah
Salt Lake City, Utah 84112

Response to Reviewers

Reviewer #1

These protocols are for building a toy self-driving lab, based on LED blinking lights and a light sensor. The protocols cover hardware and software. The result in Fig. 17 is great! I think this demo is fantastic for teaching students about self-driving labs and Bayesian optimization, and its advantage over random search or classical experimental design methods (for example, the grid search you have here). As well as for educating the public, as the concept here is easy to understand. So, the protocols are certainly deserving of the pages of STAR Protocols.

Thank you for your positive feedback! We are glad you like the results in Fig. 17, and we agree that it can be useful both in classroom settings and for the public.

This is not easy for me to assess without actually having the hardware and going through the motions. But, I generally found the description of the protocols to be clear, and I think I could reproduce them. The hardware is listed and easy to purchase. Below are a few minor comments to possibly help the authors improve their protocols article.

- Why do we need to control from the cloud? Can't we run this on our local computers and avoid the extra steps of hooking it up to the cloud?

A note was added to the Additional Prerequisites section to clarify this:

The purpose of this rather than using a hardwired connection is to emphasize the notion of "cloud experimentation", where the host and the client may be separated by large geographical distances. For more context, see <https://github.com/sparks-baird/self-driving-lab-demo/discussions/91> and <https://github.com/sparks-baird/self-driving-lab-demo/discussions/62>. For links to a simple example using a wired connection, see Problem 2:

The commentary in the Troubleshooting section has also been updated:

Problem 2: Can I use this without connecting to the internet?

Potential solution:

A simple example of wired communication between a computer and the microcontroller for the microcontroller host code and a Jupyter notebook tutorial (client) can be found at <https://github.com/sparks-baird/self-driving-lab-demo/tree/main/src/extra/nonwireless> [permalink] and <https://github.com/sparks-baird/self-driving-lab-demo/blob/main/notebooks/5.0-nonwireless-search.ipynb> [permalink], respectively. While possible with some modification, data communication via a USB cable is not actively supported for new releases of microcontroller host code nor the advanced tutorials. For private, secure, wireless communication between the Pico W microcontroller and the client (e.g., Jupyter notebook running locally), a free, private HiveMQ instance can be set up per the instructions in Software Setup. For recommendations regarding connecting to a 2.4 GHz network (e.g., in university classroom settings) see <https://github.com/sparks-baird/self-driving-lab->

demo/discussions/83 and <https://github.com/sparks-baird/self-driving-lab-demo/discussions/88>. See also page 3.

- Clarify up front, the sculpting wire is just a way to position the sensor, and will not be used to conduct electricity (correct?).

This has been clarified.

- Make the YouTube video more prominent in the article, as I think this is much more natural than a PDF (no offense to STAR Protocols)...

This is a great suggestion. It has been added to both the summary and the data availability sections.

- In the key resource table, instead of or in addition to the abbreviations, can you write out what the piece is, in plain terms? And why it is needed/its purpose? For example, CONN HEADER VERT. CBL USB2.0. are opaque names.

This is a great point. The abbreviations have largely been removed and replaced with more readable entries. Since the latest version includes some ANDs and ORs, to avoid extra confusion, the purpose has been left out. I'd prefer to add another column describing the purpose, but I think it needs to follow the key resource table format.

- The images for example in steps 8-10 are very grainy and poor quality.

The image resolution has been increased where feasible.

- For Bayesian optimization, is the code for this included?

A link has been added to both the analysis section and the data and code availability section.

Reviewer #2

IN STAR-PROTOCOLS-D-23-00017 Baird and Sparks provide a wonderful example of a low cost tool that will allow educators to introduce their students to the concepts of a self-driving lab. The manuscript is lovingly written (and the accompanying YouTube video is excellent even if it does illustrate that the authors have relatively poor musical sensibilities).

Thank you for the positive feedback! We are glad you enjoyed the manuscript and the video, and the music stems in part from Sterling's past life as a breakdancer. Thanks for bearing through it 😊

The authors carefully followed the Protocol Template with section name and timings. They included a Key Resources table that included direct links to a digikey order. I can't guarantee that those links work forever, but for right now they function.

Agreed, we hope the links will last, and we will try to check back periodically to ensure working links.

The steps (as written) are reasonably easy to follow. I personally find that words don't do a great job of describing physical builds and that their manuscript could have been more like Ikea instructions to promote clarity. I would suggest that the YouTube video be explicitly called out as I found that to be very instructive.

I think making the build instructions more like Ikea instructions is a great idea. We will try to incorporate this into future manuscripts. Great suggestion about including the video. We have added that as the last sentence to the Summary section to make it prominent as well as to the data and code availability section.

I really only have two major concerns for this protocol:

1. Folks will have versions of Python on their computer already and there may be some hesitancy around downloading Thonny and a new version of Python. I tried pip installing thonny and it did some things with my packages which may or may not have broken my Python Environment. Is there a MicroPython IDE that is compatible with Anaconda? This will make folks like me more comfortable. This might be a good split point (or an opportunity to remind folks that they should have (and use) an experimental environment for things like this.

This is a great point. From what I can tell, when Thonny is installed, it installs its own Python version (for mine, it is installed at `C:\Users\<username>\AppData\Local\Programs\Thonny\python.exe`.

I use Miniconda instead of Anaconda distribution, so the package conflicts may not appear on mine, but I gave the `pip` installation procedure a try, and it seems to work OK in a fresh conda environment:

```
conda create -n sdl-demo-thonny python==3.10.*  
conda activate sdl-demo-thonny  
pip install self-driving-lab-demo  
thonny
```

Additional content has been added related to this.

2. I found sections on MongoDB and HiveMQ confusing in that they are labeled optional but the troubleshooting makes them almost seem mandatory? Likewise the YouTube video seems to make it mandatory. This part could be further clarified.

Thanks for the great suggestion. This has been clarified in both places as follows:

- a. (Optional) Set up a MongoDB database backend. If ignored, the demo will function, just without logging data to a database (i.e., the user becomes responsible for saving the data on the client side).
...
b. (Optional) Create your own HiveMQ instance. If this setup is ignored, the demo will function properly; however, the hardware commands and sensor data will be transmitted via a default HiveMQ instance for which the credentials are public. Setting up your own HiveMQ instance ensures that the data you transfer remains private and secure. Other MQTT brokers such as Mosquitto or Adafruit IO are available. At the time of writing, we recommend HiveMQ because it provides free instances with generous limits. Setting up a private MQTT broker is in line with best practices for internet of things (IoT) security.

At the time the video was made, the MongoDB data logging was not implemented, and `test.mosquitto.org` was being used as the MQTT broker, meaning anyone could be listening in and even sending commands to

the device. The switch to HiveMQ was to allow for a free way to set up a private MQTT broker in order to follow [best practices for IoT security](#).

I think this is a really awesome project and it is unfortunate that I couldn't get all of the pieces together in time to have my students do an en masse build. I will continue to interact with the authors moving forward and can hopefully provide them

Thank you for the great review, and we look forward to our future interactions as well!



Click here to access/download
Graphical Abstract
graphical-abstract.jpg

[Click here to view linked References](#)

Building a “Hello World” for Self-driving Labs: The Closed-loop Spectroscopy Lab Light-mixing Demo (CLSLab:Light)

Sterling G. Baird^{1,2,*} and Taylor D. Sparks^{1,3,**}

¹Materials Science & Engineering Department, University of Utah, Salt Lake City UT USA, 84108

²Technical contact

³Lead contact

*Correspondence: sterling.baird@utah.edu

**Correspondence: sparks@eng.utah.edu

Summary

Learn how to build a Closed-loop Spectroscopy Lab: Light-mixing Demo (CLSLab:Light) to perform color-matching via RGB LEDs and a light sensor for under 100 USD and less than an hour of setup. Our tutorial covers ordering parts, verifying prerequisites, software setup, sensor mounting, testing, and an optimization algorithm comparison tutorial. We use secure IoT-style communication via MQTT, MicroPython firmware on a pre-soldered Pico W microcontroller, and the self-driving-lab-demo Python package. A video tutorial is available at <https://youtu.be/D54yfxRSY6s>.

For complete details on the use and execution of this protocol, please refer to Baird et al.⁽¹⁾.

Graphical abstract

Before you begin

The protocol below describes how to set up Closed-loop Spectroscopy Lab: Light-mixing Demo (CLSLab:Light), a “Hello, World!” for a “self-driving” (i.e. autonomous) laboratory (SDL)⁽²⁾ using a Pico W microcontroller, LEDs, a light sensor, and Bayesian optimization. CLSLab:Light incorporates key principles for SDLs including sending commands, receiving sensor data, physics-based simulation, and advanced optimization. This “Hello, World!” introduction is accessible to students, educators, hobbyists, and researchers for less than 100 USD, a small footprint, and under an hour of setup time. For a full video build tutorial, please refer to <https://youtu.be/D54yfxRSY6s>. There are some deviations between the instructions in the YouTube video build tutorial and recent versions of the self-driving-lab-demo Python package. In particular, see steps 13 and 14.

Order Required Parts

Timing: 5 min (not including shipping time)

1. Order the parts: (<https://www.digikey.com/short/qztj2jt7> AND [Pico W with pre-soldered headers](#)) OR <https://www.digikey.com/short/vtzjbvr2>. A visual summary of parts is given in Figure 1.

Note: For the first option, the total is 68.61 USD (or 73.72 USD including optional parts) + shipping as of 2022-03-06.

Note: The authors plan to periodically check and update the “DigiKey Order” link at <https://hackaday.io/project/186289-autonomous-research-laboratories> in case of part shortages or deprecation.

Note: In case of part shortages, many products may also be found on the Adafruit website.

Critical: If you’d like to avoid soldering, you will need to source a Pico W with headers or a Pico WH separately, such as PiShop’s [Pico W’s with pre-soldered headers](#). See also [Raspberry Pi’s supported resellers for the Pico W](#).

Note: The sculpting wire needs to be 14 gauge (2 mm) or thinner, including the insulation jacket, and rigid enough to support the sensor. The sculpting wire is only used for mounting purposes, not to conduct electricity. Sculpting wire is [also available at Amazon](#). Approximately 3' is required. See [Problem 5](#):

Note: The purpose of the wall adapter is so that, after initial setup, the demo can be powered standalone where communication happens purely via Wi-Fi.

Note: The hardware and software was designed to work with the Pico W, though the setup can be adapted for other microcontrollers. See [Problem 1](#):

Note: The bill of materials, not including the sculpting wire, is also [available at Adafruit](#).

Figure 1

Additional Prerequisites

Timing: N/A

2. **Critical:** Ensure access to a 2.4 GHz Wi-Fi network (SSID + password)

Note: The purpose of using a wireless connection rather than a hardwired one is to capture

the principles behind “cloud experimentation”, where the host and the client may be separated by large geographical distances. Additionally, this allows for a computer to only be required for initial setup such that the device can function standalone, waiting to receive commands and send sensor data. This captures best practices of a scaled-up cloud-accessible lab or network of labs. For more context, see <https://github.com/sparks-baird/self-driving-lab-demo/discussions/91> and <https://github.com/sparks-baird/self-driving-lab-demo/discussions/62>. For links to a simple example using a wired connection and related discussion, see **Problem 2:**

Note: The Pico W only supports 2.4 GHz Wi-Fi networks. See [self-driving-lab-demo #76](#) for additional context and recommendations on setting up a 2.4 GHz Wi-Fi network, if not already available.

Note: WPA enterprise networks such as Eduroam and other networks that use captive portals (most schools, coffee shops, etc.) are not yet supported by MicroPython. It needs to be a network such that on a computer, you can click on the Wi-Fi name (SSID), enter the password, and click connect (no additional steps). Check to see if your institution offers network support for internet of things devices (e.g., [ULink at University of Utah](#)).

Note: Home networks can have both a 5G and a 2.4 GHz network (e.g. “My Network 5G” and “My Network”)

Critical: If you use a mobile hotspot, you may need to use your device’s “extended compatibility” feature to drop the mobile hotspot from 5G to 2.4 GHz. See also [prepaid, long-expiry hotspot](#) and [classroom demos with standalone network access](#) discussions, which includes a summary of recommendations for prepaid mobile hotspots

3. Ensure access to a computer (for initial setup only)

Note: At a minimum, the computer needs to be able to run the Thonny editor (lightweight) and it must have at least one USB-A port

4. If the headers are not already soldered onto the microcontroller, ensure access to a soldering iron and soldering wire (thinner is better in this case)
5. **Optional:** Ensure the Pico W can successfully connect to a computer
 - a. Hold the BOOTSEL button on the Pico W while connecting the Pico W to your computer via the USB cable.

Note: If a new drive appears, that indicates that the Pico W is working normally

Note: If soldering, be careful only to heat the gold pads to avoid damaging the circuitry

Key resources table

| REAGENT or RESOURCE | SOURCE | IDENTIFIER |
|---|---|---|
| Deposited Data | | |
| Red, Green, and Blue LED Spectral Data | Baird, S. G.; Sparks, T. D. What Is a Minimal Working Example for a Self-Driving Laboratory? Matter 2022, 5 (12), 4170–4178. https://doi.org/10.1016/j.matt.2022.11.007 . | https://github.com/sparks-baird/self-driving-lab-demo/tree/v0.8.2/src/self_driving_lab_demo/data |
| Software and Algorithms | | |
| self-driving-lab-demo v0.8.2 | Baird, S. G.; Sparks, T. D. What Is a Minimal Working Example for a Self-Driving Laboratory? Matter 2022, 5 (12), 4170–4178. https://doi.org/10.1016/j.matt.2022.11.007 . | https://github.com/sparks-baird/self-driving-lab-demo |
| Other | | |
| AS7341 Color Sensor | DigiKey (Adafruit Product) | Cat#1528-4698-ND |
| Grove to Stemma-QT adapter | DigiKey (Adafruit Product) | Cat#1528-4528-ND OR Cat#1528-4528-ND |
| Raspberry Pi Pico W with pre-soldered headers OR (Raspberry Pi Pico W AND Header pins with 20 positions and 2.54 mm pitch (x2)) | PiShop OR (DigiKey-Adafruit Product AND DigiKey-Amphenol CS) | Cat#ASM-1918 OR (Cat#2648-SC0918CT-ND AND Cat#10129378-920001BLF-ND) |
| USB-A to USB-B Cable | DigiKey (Adafruit Product) | Cat#380-1431-ND |
| Maker Pi Pico base (without Pico) | DigiKey (Adafruit Product) | Cat#3614-MAKER-PI-PICO-NB-ND |
| AC/DC Wall Mount Adapter 5V 5W | DigiKey (Adafruit Product) | Cat#1470-2768-ND |
| 18 AWG Hook-up solid black wire, 100' (Outer diameter 14 AWG or higher) | DigiKey (Remington Industries) | Cat#2328-18UL1007SLDBLA-ND |
| Terminal Binding Post M2.5 (Optional) | DigiKey (Keystone Electronics) | Cat#36-8737-ND |
| 128MB Micro SD Memory Card (Optional) | DigiKey (Adafruit Product) | Cat#1528-5250-ND |

Step-by-step method details

Hardware Setup

Timing: 20 min

Unless pre-soldered, attach the headers onto the Pico W, mount the light sensor so that the pinhole is facing the red green blue (RGB) LED, connect the light sensor to the board, and get the microcontroller ready for firmware installation.

1. Unless pre-soldered, Solder headers onto the Pico W or [use a hammer header pin install rig for Pico W](#)

Note: If soldering, insert the Pico W headers into the Maker Pi Pico base, place the Pico W on top of the headers, and solder the headers to the Pico W ([MagPi guide](#), [Tom's hardware guide](#), or [YouTube video](#)), and remove the Pico W from the Maker Pi Pico base

Note: Pico install rigs are not compatible with the Pico W. It must be labeled explicitly as "Pico W".

2. Prepare 3 feet of sculpting wire (cut with wire cutters or bend until it breaks)
3. Thread the sculpting wire through each mounting hole on the Maker Pi Pico base, then twist the wires together near the RGB LED. See Figure 2, Figure 3, and Methods Video S1.

Note: This setup will allow the position and orientation of the sensor to be both adjustable and steady.

Figure 2

Figure 3

4. Continue twisting until you have 4 to 6 inches of twisted wire, and ensure that there are at least 3 inches of loose, untwisted wire at each end. See Figure 2 and Methods Video S1.

Note: (the leftover, untwisted wire will be threaded through the mounting holes of the light sensor in the next step). For a more modular alternative of fixturing the wire ends to the Maker Pi Pico base, see [Problem 5](#):

5. Thread the same sculpting wire through the AS7341 light sensor and position the sensor so the pinhole is facing approximately 3 to 4 inches away from the RGB LED. See Figure 4 and Methods Video S2.

Figure 4

6. Connect the Grove/Stemma-QT connector into Grove port 6 (GP26&27) and the AS7341, insert the SD card (**optional**), insert the Pico W, and while holding the BOOTSEL button, connect the Pico W to the computer. See Figure 5 and Methods Video S3.

Figure 5

Software Setup

Timing: 20 min

Install the MicroPython firmware onto the Pico W microcontroller, enter the Wi-Fi credentials, and upload the source code files.

7. Download and install [Thonny](#), a Python IDE with native support for microcontrollers, onto your computer. See Methods Video S4.
 - a. Choose the platform appropriate for you (in my case, this is Windows 64-bit, Python 3.10).
 - b. When installing, use the default settings: "Standard (default)". Thonny comes with its own version of Python located by default at C:\Users\<username>\AppData\Local\Programs\Thonny\python.exe on Windows computers.
 - c. It is not anticipated that this will cause conflicts with existing installations of Python; however, for conda users, an isolated installation may be performed via the following commands in a conda shell:
`conda create -n sdl-demo-thonny python==3.10.*`
`conda activate sdl-demo-thonny`
`pip install thonny`
`thonny`
8. Click on the lower-right dropdown and click "Install MicroPython", which will install the microcontroller firmware onto the Pico W. See Figure 6 and Methods Video S4.

Figure 6

9. Choose "MicroPython variant: Raspberry Pi - Pico W / Pico WH" and click install. See Figure 7 and Methods Video S4.

Figure 7

10. Change the interpreter from Local Python 3 to MicroPython (Raspberry Pi Pico), which will open a shell that can be used to enter MicroPython commands that run directly on the Pico W. See Figure 8 and Methods Video S4.

Figure 8

11. In Thonny's menubar, click "View" then "Files" to open a sidebar which shows both your local computer's files (top) and the files on the Pico W (bottom). See Figure 9 and Methods Video S5.

Figure 9

12. Download *sdl_demo.zip* from [the latest release at self-driving-lab-demo](https://self-driving-lab-demo) to your computer and unzip it. See Methods Video S5.
13. In Thonny, navigate to the unzipped *sdl_demo* folder, open *secrets.py*, enter your Wi-Fi network name (SSID) and password as Python strings, and save *secrets.py*. See Figure 10, Figure 11, and Methods Video S5.

Optional: you can create your own MongoDB Atlas database and enter values for MONGODB_API_KEY, MONGODB_COLLECTION_NAME, and DEVICE_NICKNAME into *secrets.py* (see below).

Optional: you can create your own HiveMQ instance and enter *secrets.py* credentials for HIVEMQ_USERNAME, HIVEMQ_PASSWORD, and HIVEMQ_HOST (see below).

Figure 10

Figure 11

- a. Set up a MongoDB database backend.

Note: If ignored, the demo will function, just without logging data to a database (i.e., the user becomes responsible for saving the data on the client side). See **Problem 3:**

- i. Create an account at
<https://www.mongodb.com/cloud/atlas/register>
- ii. Create a free, Shared Cluster. See Figure 12.

Note: optionally rename Cluster0 to something of your choice, e.g. self-driving-labs. You can leave the default provider as-is.

Figure 12

- iii. Navigate to “Data Services” → “Deployment” → “Database” and click “Browse Collections”. See Figure 13.

Figure 13

- iv. Click “Add My Own Data”
- v. Enter a database name (e.g., clslab-light-mixing) and collection name (e.g., test).
- vi. Copy the names into MONGODB_DATABASE_NAME and MONGODB_COLLECTION_NAME in *secrets.py*.
- vii. Navigate to “Data Services” → “Services” → “Data API”, use the dropdown to select your cluster, and click “Enable Data Access from

the Data API". See Figure 14.

Figure 14

- viii. Note the app name in the “URL Endpoint” box of the form “<https://data.mongodb-api.com/app/<data-abc123>/endpoint/data/v1>” where <data-abc123> is the app name. See Figure 15.

Figure 15

- ix. Copy the app name into the MONGODB_APP_NAME variable in *secrets.py*.
- x. Click “Create API Key”, enter a name of your choice (e.g. clslab-light), and click “Generate API key”. See Figure 16.

Figure 16

- xi. Copy the API key and store it somewhere secure, then paste the API key into the MONGODB_API_KEY variable in *secrets.py*.
- b. Create your own HiveMQ instance.

Note: If this setup is ignored, the demo will function properly; however, the hardware commands and sensor data will be transmitted via a default HiveMQ instance for which [the credentials are public](#). Setting up your own HiveMQ instance ensures that the data you transfer remains private and secure. Other MQTT brokers such as [Mosquitto](#) or [Adafruit IO](#) are available. At the time of writing, we recommend HiveMQ because it provides free instances with generous limits. Setting up a private MQTT broker is in line with [best practices for internet of things \(IoT\) security](#) and should be used especially when working with sensitive data.

- i. Navigate to <https://www.hivemq.com/mqtt-cloud-broker/>, click “Try out for free”, and create an account
- ii. Set up credentials by entering a username and password and press “ADD”. See Figure 17.

Figure 17

- iii. Navigate to the “Clusters” tab and copy the URL (e.g., abc123.s2.eu.hivemq.cloud) to HIVEMQ_HOST in *secrets.py*. Also update HIVEMQ_USERNAME and HIVEMQ_PASSWORD with the

username and password from the previous step. See Figure 18.

Figure 18

- iv. Create a certificate using the Google Colab notebook at <https://github.com/sparks-baird/self-driving-lab-demo/blob/v0.7.3/notebooks/7.2.1-hivemq-openssl-certificate.ipynb>
Note: This file is used to do secure authentication via HiveMQ.
- v. Enter the server address (i.e., HIVEMQ_HOST) into *secrets.py* and run the Google Colab cells
- vi. Follow the instructions to download the hivemq-com-chain.der file to the unzipped sdl_demo folder.

14. Upload files to the Pico W microcontroller. See Figure 19 and Methods Video S6.

- a. While holding Ctrl (Windows) or Cmd (Mac), select "lib", "main.py", "hivemq-com-chain.der", and "secrets.py"

Note: hivemq-com-chain.der is not mentioned in the YouTube tutorial, as it was not implemented at the time of creating the video.

- b. Right click in the gray region
- c. Click "Upload to /"

Figure 19

15. Double click to open *main.py*, click the green play button (i.e., run the code on the Pico W), and note the PICO ID that prints to the command window ("prefix/picow/<PICO_ID>/"). See Figure 20 and See Methods Video S6.

Critical: *main.py* needs to be run on the Pico W microcontroller (host), not on your local machine (client).

Note: This will act as the “password” to control the demo.

Figure 20

Control from the cloud

Timing: 10 min

Control the device via internet-of-things style communication (MQTT) and run a basic optimization comparison of grid search vs. random search vs. Bayesian optimization.

16. [Open notebooks/4.2-paho-mqtt-colab-sdl-demo-test.ipynb in Google Colab](#). See Methods Video S7.

17. Scroll to the first code cell and click the play button to install the self-driving-lab-demo Python package. See Figure 21 and Methods Video S7.

Figure 21

18. Copy the PICO ID from the Thonny editor and paste it in place of "test" (without quotes). See Figure 22. An example image of the output is given in Figure 23. See also Methods Video S8.

Note: the actual output to the command window may vary in future releases.

Note: If you leave PICO_ID set to "test", this will control a public demo maintained by the authors for testing and demonstration purposes. The authors will strive to keep this public test demo available for the foreseeable future with minimal downtime.

Figure 22

Figure 23

19. Run the remaining code cells. See Methods Video S8, Methods Video S9, and Methods Video S10.
 - a. Instantiate a SelfDrivingLabDemo class
 - b. Perform optimizations for grid search, random search, and Bayesian optimization

Expected outcomes

It is expected that users will successfully set up the hardware and software for a closed-loop experiment. Further, users will run their first "autonomous drive" given in an example interactive notebook and explore [additional example notebooks](#).

Figure 24 shows a comparison of optimization results for grid search vs. random search vs. Bayesian optimization averaged over repeat campaigns with standard deviation error bands, where Bayesian optimization, on average, performs the best. Figure 25 shows one of the outputs from the cloud-based control notebook of best error so far vs. iteration number comparing grid search vs. random search vs. Bayesian optimization. Typically, grid search is the least efficient, Bayesian optimization is the most efficient, and random search is somewhere in-between. Figure 26, Figure 27, and Figure 28 show the points that were searched for a given campaign for grid search, random search, and Bayesian optimization, respectively. Finally, Figure 29 shows the true, underlying target color (defined by red, green, and blue values) and the best parameter set based on minimizing error between the observed spectrum and the target spectrum for each of the optimization methods.

Figure 24

Figure 25

Figure 26

Figure 27

Figure 28

Figure 29

Quantification and statistical analysis

Discrete Fréchet distance, as implemented in https://github.com/cjekel/similarity_measures, is used to assess the mismatch between the currently observed spectrum and the target spectrum, where the target spectrum is determined by arbitrarily choosing a random set of RGB values and measuring the sensor data for the fixed, random set of RGB values. Lower Fréchet distances correspond to better matches between the observed and target spectra (i.e. lower error).

An example JSON document logged to a MongoDB database backend containing experimental data for a single run is given as follows:

```
{
    "utc_timestamp": "2022-11-4 06:51:16",
    "ch510": 354,
    "ch620": 5671,
    "ch410": 188,
    "ch440": 3675,
    "ch583": 2756,
    "_input_message": {
        "_session_id": "542e6e80-9c50-4c41-95a5-832603b96238",
        "B": 31,
        "atime": 100,
        "gain": 128,
        "astep": 999,
        "_experiment_id": "9b50c819-db8f-476f-b601-dbe79e871a46",
        "G": 3,
        "integration_time": 280.78,
        "R": 41,
    },
    "onboard_temperature_K": 294.1085,
    "sd_card_ready": True,
    "ch470": 2827,
    "ch550": 498,
    "ch670": 277,
}
}
```

The experimental parameters for two JSON documents are given in Table 1.

Table 1. Example of data obtained from two experiments. The LED parameters are red (R), green (G), blue (B). The sensor settings are atime, gain, astep (affects integration time and intensity). The measured output values are of the form "ch###" where the three digit number corresponds to the full-width half-max (FWHM) wavelength being measured.

| utc_timestamp | onboard_temperature_K | R | G | B | atime | gain | astep | ch410 | ch440 | ch470 | ch510 | ch550 | ch583 | ch620 | ch670 |
|----------------|-----------------------|----|---|----|-------|------|-------|-------|-------|-------|-------|-------|-------|-------|-------|
| 11/4/2022 6:40 | 292.7041 | 41 | 3 | 31 | 100 | 128 | 999 | 188 | 3674 | 2828 | 354 | 498 | 2748 | 5661 | 276 |
| 11/4/2022 6:51 | 294.1085 | 41 | 3 | 31 | 100 | 128 | 999 | 188 | 3675 | 2827 | 354 | 498 | 2756 | 5671 | 277 |

The code for grid search, random search, and Bayesian optimization is hosted at https://github.com/sparks-baird/self-driving-lab-demo/blob/main/src/self_driving_lab_demo/utils/search.py [permalink].

Limitations

Environmental noise (e.g. light conditions) and hardware variation (LED, sensor, sensor positioning, etc.) may affect the results obtained.

Troubleshooting

See the [GitHub issue tracker](#) for existing known issues or to post a new issue. See the [GitHub discussions](#) for general questions and discussion.

Problem 1:

Can I use this with alternate microcontrollers or firmware?

Potential solution:

The hardware configuration and software were designed based on Raspberry Pi's Pico Wireless (Pico W) microcontroller. Libraries exist for LED control and the AS7341 light sensor in CircuitPython and Arduino. The hardware and configuration and software can be adapted for other microcontrollers.

Contributions at <https://github.com/sparks-baird/self-driving-lab-demo/> are welcome. See [Order Required Parts](#).

Problem 2:

Can I use this without connecting to the internet?

Potential solution:

A simple example of wired communication between a computer and the microcontroller for the microcontroller host code and a Jupyter notebook tutorial (client) can be found at <https://github.com/sparks-baird/self-driving-lab-demo/tree/main/src/extra/nonwireless> [permalink] and <https://github.com/sparks-baird/self-driving-lab-demo/blob/main/notebooks/5.0-nonwireless-search.ipynb> [permalink], respectively. While possible with some modification, data communication via a USB cable is not actively supported for [new releases of microcontroller host code](#) nor [the advanced tutorials](#). The status of this feature is being tracked at <https://github.com/sparks-baird/self-driving-lab-demo/issues/193>. For private, secure, wireless communication between the Pico W microcontroller and the client (e.g., Jupyter notebook running locally), a free, private HiveMQ instance can be set up per the instructions in [Software Setup](#). For recommendations regarding connecting to a 2.4 GHz network (e.g., in university classroom settings) see <https://github.com/sparks-baird/self-driving-lab-demo/discussions/83> and <https://github.com/sparks-baird/self-driving-lab-demo/discussions/88>. See also step [Additional Prerequisites](#).

Problem 3:

Can I use this without logging to a MongoDB backend?

Potential solution:

If the MongoDB credentials are left to their default dummy values in secrets.py, then logging to the MongoDB backend will fail and the device will simply notify the user rather than exit the program. In other words, the device will function normally without database logging. The same applies for logging to an onboard SD card. If an SD card is detected, the microcontroller will write backup data to it, otherwise this step will be skipped. See step 13.a.

Problem 4:

The Stemma-QT to Grove connectors (or other items) are out-of-stock.

Potential solution:

First, look at Adafruit and other vendors to see if it is available. Note that Cat#[1528-4424-ND](#) is incompatible with the Maker Pi Pico base due to the adapter housing blocking it from being plugged in fully. If no Stemma-QT to Grove connectors can be located, another alternative is using a Stemma-QT to header pin cable (DigiKey Cat#[1528-4209-ND](#)) and plugging directly into the GPIO pins that correspond to Grove Port #6 of the Maker Pi Pico base. For other items that may be out of stock on DigiKey or Adafruit, other vendors may be used (e.g., [AS7341 light sensor from electromaker](#)). See

Order Required Parts.

Problem 5:

The sculpting wire doesn't fit through the mounting holes.

Potential solution:

Ensure that the outer diameter of the sculpting wire is 14 AWG or higher (i.e., 1.628 mm or thinner). Enameled wire (often advertised as sculpting wire) has a very thin coating, whereas electrical wiring typically has a non-negligible insulation thickness. Optionally, for a more modular setup, a single M2.5 binding post ([Digikey Cat#36-8737-ND](#)) can be used to clamp the wire via a single mounting hole instead of looping the wire through each of the mounting holes. See steps [Order Required Parts](#) and step 3.

Problem 6:

My SD card isn't being recognized.

Potential solution:

First, we note that use of the micro SD card is optional and serves the purpose of onboard backup data logging. The 128 MB micro SD card recommended in this work (DigiKey [Cat#1528-5250-ND](#)) has been tested with the rest of the components. First, try removing the micro SD card completely and reinsert it, making sure there is an audible "click". If the microcontroller fails to detect the micro SD card, then there may be a defect in the micro SD card or the Maker Pi Pico base. Try ordering an extra micro SD card (same one recommended above), and if it suddenly works, you should be able to request a refund on the first SD card. If it still does not work, contact the seller of the Maker Pi Pico base to request a replacement. If not using the recommended SD card, the card formatting may be incompatible with MicroPython (see <https://github.com/CytronTechnologies/MAKER-PI->

[PICO/issues/4](#)). In this case, you will likely need to purchase a different type of SD card. See [Order Required Parts](#) and step 3.

Resource availability

Lead contact

Further information and requests for resources and reagents should be directed to and will be fulfilled by the lead contact, Taylor D. Sparks sparks@eng.utah.edu.

Materials availability

This study did not generate new unique reagents.

Data and code availability

The datasets and code generated during this study are available on GitHub:

<https://github.com/sparks-baird/self-driving-lab-demo>. The recommended option for ordering parts is <https://www.digikey.com/short/qztj2jt7> AND [Pico W with pre-soldered headers](#) to avoid soldering. Alternatively, <https://www.digikey.com/short/vtzjbvr2> is a standalone DigiKey order, but requires soldering headers onto the Pico W. For a full video build tutorial, please refer to <https://youtu.be/D54yfxRSY6s>. Code for grid search, random search, and Bayesian optimization is hosted at https://github.com/sparks-baird/self-driving-lab-demo/blob/main/src/self_driving_lab_demo/utils/search.py. A version of record for v0.8.2 is given in DOI: [10.5281/zenodo.7855493](https://zenodo.10.5281/zenodo.7855493). To cite all versions, see DOI: [10.5281/zenodo.7855492](https://zenodo.10.5281/zenodo.7855492). We encourage readers to use the latest version of the self-driving-lab-demo Python package and revert to v0.8.2 if breaking changes occur that prevent use of the latest package with the instructions in this protocol. We will make efforts to minimize changes that are backwards incompatible, and we would highly appreciate if users would check [the existing issues](#) (see both open and closed issues) and [open a new issue](#) in the [GitHub issue tracker](#) if not already present in the existing issues. A free GitHub account can be created to comment on existing issues or open new issues, and a [GitHub markdown syntax guide](#) is available. Alternatively, users may also contact the authors at the emails listed in the author affiliations.

Acknowledgments

This work was supported by the National Science Foundation under Grant No. DMR-1651668.

Author contributions

Sterling G. Baird: Conceptualization, Methodology, Software, Writing – Original Draft, Writing – Review & Editing, Visualization, Taylor D. Sparks: Supervision, Funding Acquisition

Declaration of interests

The authors have been exploring selling at-cost kits via a crowdfunding platform called [GroupGets](#). As of 2023-04-23, the cost-breakdown is as follows: List price: 80 USD + shipping (depends on location). GroupGets fee: 15.57 USD (10%+2.9%+5 USD). Hardware cost: 51.19 USD. Labor/testing: 20 min (est. 6 USD). Profit: 2.24 USD, offset somewhat by bulk pricing discounts. Note that this assumes purchasing a Pico W with presoldered headers. See [round 1](#) and [round 2](#) and a discussion of [packaging open-source hardware as commercial kits](#).

References

- a. S.G. Baird, T.D. Sparks, “What is a Minimal Working Example for a Self-driving Laboratory?” Matter, Cell Press, 2022. 5 (12), 4170–4178. <https://doi.org/10.1016/j.matt.2022.11.007>.
- b. Seifrid, M.; Hattrick-Simpers, J.; Aspuru-Guzik, A.; Kalil, T.; Cranford, S. Reaching Critical MASS: Crowdsourcing Designs for the next Generation of Materials Acceleration Platforms. Matter, Cell Press, 2022. 5 (7), 1972–1976. <https://doi.org/10.1016/j.matt.2022.05.035>.

Figure legends

- Figure 1: Visual bill of materials
- Figure 2: Wire mounting instructions
- Figure 3: Wire mounting schematic
- Figure 4: Light sensor mounting instructions
- Figure 5: Hardware connections
- Figure 6: Firmware installation dropdown
- Figure 7: MicroPython installation dialogue box
- Figure 8: Interpreter dropdown
- Figure 9: Opening the files sidebar
- Figure 10: Editing secrets.py
- Figure 11: Saving secrets.py
- Figure 12: Setting up a MongoDB shared cluster
- Figure 13: Create a MongoDB database
- Figure 14: Enable the Data API
- Figure 15: Retrieve MONGO
- Figure 16: Create Data API key
- Figure 17: Set up HiveMQ credentials
- Figure 18: Locate the HiveMQ host URI
- Figure 19: Uploading source files to microcontroller
- Figure 20: Running main.py
- Figure 21: Python package installation
- Figure 22: Copying the Pico ID from the Thonny editor
- Figure 23: Pasting the Pico ID into the Google Colab form box
- Figure 24: Example optimization comparison between grid search, random search, and Bayesian optimization averaged over repeated campaigns. Lower Fréchet distance between observed and target spectra is better.

Figure 25: Example optimization comparison between grid search, random search, and Bayesian optimization. Lower error is better.

Figure 26: Twenty-seven grid search points colored by the Fréchet distance between the target spectrum and the sensor data evaluated at each grid point.

Figure 27: Twenty-seven random search points colored by the Fréchet distance between the target spectrum and the sensor data evaluated at each grid point.

Figure 28: Twenty-seven Bayesian optimization points colored by the Fréchet distance between the target spectrum and the sensor data evaluated at each grid point.

Figure 29: The true, underlying RGB target (purple diamond) and the best observed points for grid search (blue circle), random search (red circle), and Bayesian optimization (green circle). Bayesian optimization gave the closest match to the true target.

Methods Video S1: Thread the mounting wire through the mounting holes of the Maker Pi Pico base. See step 3 and step 4.

Methods Video S2: Thread the remaining mounting wire through the mounting holes of the AS7341 light sensor and position the sensor above the LEDs. See step 5.

Methods Video S3: Attach the Pico W and the AS7341 light sensor to the Maker Pi Pico base, then connect the USB cable from the Pico W to the computer while holding down the BOOTSEL button. See step 6

Methods Video S4: Download the Thonny editor and install the MicroPython firmware onto the Pico W. See steps 7, 8, 9, and 10.

Methods Video S5: Download the source code from GitHub, unzip it, and enter Wi-Fi credentials. See steps 11, 12, and 13.


Methods Video S6: Upload the source code to the Pico W and run the main.py script. See steps 14 and 15.


Methods Video S7: Open the cloud-control Jupyter notebook via Google Colab and install the self-driving-lab-demo Python package. See steps 16 and 17.


Methods Video S8: Copy-paste the PICO ID from Thonny to Colab and control the setup remotely through the “evaluate” command. See steps 18 and 19.




Methods Video S9: Perform the “Hello, World!” of optimization, comparing grid search vs. random search vs. Bayesian optimization. See step 19.


Methods Video S10: Visualize the results of the optimization comparison, See step 19.











✓ Local Python 3 • Thonny's Python


[Install MicroPython...](#)

[Install CircuitPython...](#)

[Configure interpreter...](#)

Local Python 3 • Thonny's Python

Figure 7

[Click here to access/download;Figure;fig7.jpg](#)


✓ Local Python 3 • Thonny's Python

MicroPython (Raspberry Pi Pico) • COM4

MicroPython (RP2040) • COM4


Configure interpreter...

Local Python 3 • Thonny's Python



A screenshot of a code editor window titled "secrets.py". The code contains three lines of Python code. Line 1 is a multi-line string docstring. Line 2 defines a variable "SSID" with a value that is completely redacted by a large black rectangle. Line 3 defines a variable "PASSWORD" with a value that is also completely redacted by a large black rectangle.

```
1 """WiFi Login Info. Rename this file to secrets.py"""
2 SSID = "████████"
3 PASSWORD = "████████"
```



CLUSTERS > CREATE A SHARED CLUSTER

Create a Shared Cluster

Serverless Dedicated **Shared**

For learning and exploring MongoDB in a sandbox environment. Basic configuration controls. X

No credit card required to start. Upgrade to dedicated clusters for full functionality.
Explore with sample datasets. Limit of one free cluster per project.

Cloud Provider & Region AWS, N. Virginia (us-east-1) ▾

aws Google Cloud Azure

★ Recommended region ⓘ ⓘ Dedicated tier region ⓘ

NORTH AMERICA EUROPE AUSTRALIA

Figure 13

[Click here to access/download;Figure;fig13.jpg](#)

The screenshot shows the MongoDB Atlas interface for a database deployment. The left sidebar has sections for DEPLOYMENT (Database, Data Lake, PREVIEW), SERVICES (Triggers, Data API, Data Federation), and SECURITY (Quickstart, Database Access, Network Access, Advanced). The main area shows a cluster named "Cluster0" with metrics like R: 0, W: 0, Connections: 0, and Data Size: 0.0 B. A blue box highlights the "Browse Collections" button. Below the metrics is an "Enhance Your Experience" section with an "Upgrade" button. At the bottom, there's a table with columns: VERSION, REGION, CLUSTER SIZE, TYPE, BACKUPS, LINKED APP SERVICES, and ATLAS SEARCH. The VERSION is 5.0.14, REGION is AWS / N. Virginia (us-east-1), CLUSTER SIZE is M0 Sandbox (General), TYPE is Replica Set - 3 nodes, BACKUPS is Inactive, LINKED APP SERVICES is None Linked, and ATLAS SEARCH is Create Index.

Atlas Sterling's Or... Access Manager Billing All Clusters Det. Help Sterling

Data Services App Services Charts

STERLING'S ORG - 2022-03-25 > TEST

Database Deployments

Cluster0 Connect View Monitoring Browse Collections ...

Enhance Your Experience

R: 0 W: 0 Last 3d seconds Upgrade

In: 0.0 B/s Out: 0.0 B/s Last 3d seconds 0.0 B/s

Connections: 0 Data Size: 0.0 B Last 3d seconds 0.0 B

VERSION REGION CLUSTER SIZE TYPE BACKUPS LINKED APP SERVICES ATLAS SEARCH

| | | | | | | |
|--------|-------------------------------|----------------------|-----------------------|----------|-------------|--------------|
| 5.0.14 | AWS / N. Virginia (us-east-1) | M0 Sandbox (General) | Replica Set - 3 nodes | Inactive | None Linked | Create Index |
|--------|-------------------------------|----------------------|-----------------------|----------|-------------|--------------|

Figure 14

[Click here to access/download;Figure;fig14.jpg](#)

The screenshot shows the Data Services interface with the 'Data API' tab selected. The main heading is 'Enable the Data API' with the sub-instruction 'Access the data in your collection over fully managed and secure HTTPS Endpoints'. Below this, there's a dropdown menu labeled 'All Data Sources selected' and a link to 'ADVANCED CONFIGURATION'. A prominent green button at the bottom is labeled 'Enable Data Access from the Data API'. A note below the button states 'You can modify CORS settings in the configuration file'. On the left sidebar, other tabs like 'Database', 'Data Lake', and 'App Services' are visible, along with sections for 'Deployment', 'Services', and 'Security'.

test

Data Services App Services Charts

DEPLOYMENT

Database

Data Lake PREVIEW

SERVICES

Triggers

Data API

Data Federation

SECURITY

Quickstart

Database Access

Network Access

Advanced

New On Atlas

Dot

Enable the Data API

Select the data source(s) you would like to enable the API on

All Data Sources selected

ADVANCED CONFIGURATION

Enable Data Access from the Data API

You can modify CORS settings in the configuration file

Figure 15

[Click here to access/download;Figure;fig15.jpg](#)

The screenshot shows the MongoDB Data API configuration interface. The top navigation bar includes tabs for 'Data Services' (selected), 'App Services', and 'Charts'. On the left, a sidebar lists 'test' (Database), 'PREVIEW' (Data Lake), 'SERVICES' (Triggers, Data API selected), and 'SECURITY' (Quickstart, Database Access, Network Access, Advanced). The main content area is titled 'STERLING'S ORD - 2023-03-26 > TEST' and 'Data API ENABLED'. It features tabs for 'Data API' (selected), 'API Keys', 'Logs', and 'Settings'. A note states: 'By default, your Data API returns data as JSON which will cast Timestamp, ObjectId, Binary, and Decimal128 data types into strings. Visit the Settings page or add the appropriate header to change this return type to Extended JSON. Learn More'.

| URL Endpoints | Version: v1 | + | https://data.mongodb-api.com/app/[REDACTED]/endpoint/data/v1 | Copy | Advanced Settings | View our API documentation |
|---------------|------------------------------|------|---|---|-----------------------------------|--|
| Data Source | Provider/Region | Tier | Data API Access | | | |
| Cluster2 | AWS, N. Virginia (US EAST 1) | Mo | Read and Write | | | |

At the bottom, there's a 'New On Atlas' button, a 'Go to' link, and footer links for 'System Status: All Good', 'Status', 'Terms', 'Privacy', 'Atlas Blog', and 'Contact Sales'.



Create Data API Key

Create a Data API Key and be sure to store it in a secure location. You can then visit the API Key tab to [view and manage your API Keys](#) 

Configure other authentication methods for Data API in [Authentication services](#) 

Name your key

`clslab-light`

Generate API Key

-  Your Data API key only gives you access to the Data API, not direct access to data in clusters. To prevent security breaches do not distribute it to untrusted individuals or embed directly in your client applications. [Learn more about Data API keys](#) 

After you leave this page, the full private key is unavailable.

Close

Figure 17

[Click here to access/download;Figure;fig17.jpg](#)

The screenshot shows the 'Credentials' section of the HiveMQ Cloud interface. On the left, there's a sidebar with icons for 'Clusters', 'Billing', and 'Help'. The main area has tabs for 'OVERVIEW', 'ACCESS MANAGEMENT' (which is selected), 'INTEGRATIONS', 'WEB CLIENT', and 'GETTING STARTED'. A large central box is titled 'Set up credentials for your IoT devices'. It contains instructions: 'Define the credentials that your MQTT clients can use to connect to your HiveMQ Cloud cluster.' and 'Please visit the [HiveMQ documentation](#) for examples on how to use the credentials to connect an MQTT client to your cluster.' Below this, it says '(All fields are mandatory)'. There are three input fields: 'Username' (with placeholder 'At least 8 characters. Username must be unique.'), 'Password' (with placeholder 'At least 8 characters, numbers, upper- and lowercase letters.'), and 'Confirm Password' (with placeholder 'Passwords must match.'). At the bottom right of this box is a yellow 'ADD' button. To the right, a separate box is titled 'Active MQTT Credentials' with the sub-instruction 'These credentials allow MQTT clients to publish and subscribe to your HiveMQ Cloud cluster.' It contains a blue message box with a yellow exclamation mark icon and the text 'No credentials configured for your devices' followed by 'Please set up credentials to connect your IoT devices.'

Figure 18

[Click here to access/download;Figure;fig18.jpg](#)

The screenshot shows a user interface for managing clusters. On the left, there's a sidebar with icons for Clusters (selected), Billing, and Help. The main area has a title "Your Clusters" with a yellow profile icon. A "CREATE NEW CLUSTER" button is in the top right. Below it, a teal banner says "FREE Perfect for testing and small use cases" with a hexagon icon. The main content area displays a cluster entry:

| URL | PORT (TLS) |
|-------------------------------|------------|
| [REDACTED].s1.eu.hivemq.cloud | 8883 |

Below this, another row shows:

| STATUS | STARTED |
|---------|---------------------|
| Running | 12/20/2022, 2:25 PM |

A yellow "MANAGE CLUSTER" button is at the bottom right.




Figure 20

[Click here to access/download;Figure;fig20.jpg](#)

The screenshot shows the Thonny IDE interface. The title bar reads "Thonny - C:\Users\sterg\Documents\sdl_demo\main.py @ 1:55". The menu bar includes File, Edit, View, Run, Tools, and Help. The toolbar features icons for file operations and execution. The left sidebar displays a file tree under "This computer" with "C:\Users\sterg\Documents" as the root. Inside "Documents", there are several folders: Microsoft Hardware, My Kindle Content, OneNote Notebooks, overleaf, Python Scripts, and sdl_demo. The "sdl_demo" folder is expanded, showing its contents: lib, hivemq-com-chain.der, main.py, and secrets.py. The main workspace shows the Python script "main.py" with the following code:

```
1 """Run a self-driving lab on a Raspberry Pi Pico W."""
2
3 import gc
4 import json
5 import os
6 from secrets import HIVEMQ_HOST, HIVEMQ_PASSWORD, HIVE
7 from time import sleep
8
9 import ntptime
```

A "Shell" tab is visible at the bottom of the workspace.



```
▶%pip install self-driving-lab-demo
```

...

```
Requirement already satisfied: ipykernel in /usr/l
Requirement already satisfied: tqdm in /usr/local/
Requirement already satisfied: importlib-metadata
Requirement already satisfied: nbformat in /usr/lo
Requirement already satisfied: plotly in /usr/loc
Requirement already satisfied: ipywidgets in /usr/
Requirement already satisfied: numpy in /usr/local
```

The screenshot shows the Thonny IDE interface. On the left, the project tree displays files: Raspberry Pi Pico, lib, main.py (selected), and secrets.py. The main area shows Python code:

```
16 my_id = hexlify(unique_id()).decode()
17
18 prefix = f"sdl-demo/picow/{my_id}/"
19
20 print(f"prefix: {prefix}")
21
22 pixels = NeoPixel(Pin(28), 1) # one NeoPixel on Pin 28 (GP28)
```

A large callout box with a black background and white text overlays the code area, containing the instruction: "Copy the Pico ID from the Thonny editor".

The bottom section shows the terminal window with the command: `>>> %Run -c $EDITOR_CONTENT`. The output is:
prefix: sdl-demo/picow/**a123b456c789/**
Detected devices at I2C addresses: 0x20
connected.

Another callout box with a black background and white text provides an example: "For example, if you see 'prefix: sdl-demo/picow/a123b456c789/' then 'a123b456c789' (without quotes) is your Pico ID".

Figure 23

[Click here to access/download;Figure;fig23.jpg](#)

```
+ Code + Text ⌂ Copy to Drive  
from uuid import UUID  
from selfDrivingLab0 import SelfDrivingLab0Demo, mqtt_observe_sensor_data  
uniqueIdentifier = str(uuid4())  
PICO_ID = uniqueIdentifier # unique identifier  
SESSION_ID = str(uuid4()) # random session ID  
print(f"session ID: {SESSION_ID}")  
  
sd1 = SelfDrivingLab0Demo(  
    autoload=True, # perform target data experiment automatically  
    observe_sensor_data_fn=mqtt_observe_sensor_data, # (default)  
    observe_sensor_data_kwargs=dict(pico_id=PICO_ID, session_id=SESSION_ID),  
)
```

PICO_ID: [REDACTED] 35

Paste the Pico ID in place of "test" (without quotes) into the input form above and run the cell

Figure 24


[Click here to access/download;Figure;fig24.jpg](#)

Figure 25


[Click here to access/download;Figure;fig25.jpg](#)

Figure 26


[Click here to access/download;Figure;fig26.jpg](#)

Figure 27


[Click here to access/download;Figure;fig27.jpg](#)

Figure 28



[Click here to access/download;Figure;fig28.jpg](#)


Figure 29

[Click here to access/download;Figure;fig29.jpg](#)




Click here to access/download

Figure 360 and Methods Videos
1-base-threading.mp4




Click here to access/download

Figure 360 and Methods Videos
2-thread-sensor.mp4



Click here to access/download

Figure 360 and Methods Videos
3-hardware-connections.mp4




Click here to access/download

Figure 360 and Methods Videos
4-firmware-installation.mp4



Click here to access/download

Figure 360 and Methods Videos
5-source-code.mp4



Click here to access/download

Figure 360 and Methods Videos
6-source-code.mp4



Click here to access/download

Figure 360 and Methods Videos
7-cloud-control.mp4



Click here to access/download

Figure 360 and Methods Videos
8-cloud-control.mp4



Click here to access/download

Figure 360 and Methods Videos
9-optimization.mp4



Click here to access/download

Figure 360 and Methods Videos
10-visualization.mp4